



Modern Care in Motion

Supporting Your Expertise with EndoFresh® Portable Endoscopy Solution GEC-2000P



Dedicated carrier

Pack everything, including the medical display, in just 10 seconds



No tangle

3 built-in pumps with an auxiliary water-jet



No compromises

Full HD output, ample internal memory for all your needs

*GEP-2000P is currently in the process of obtaining CE approval.



Better managing infectious risks



Safeguarding vulnerable patients



Rapid response in field emergencies



Streamlining endoscopy Training



Advantages of Disposable Endoscopes



A fresh endoscope for every procedure

The sterile and single-use packaging reduces the risk of infection.



Fast setup for any clinical scenario

Always prepared to meet a wide range of medical needs.



Affordable and maintenance-free

Reliable performance without the expense or delays of repairs all your needs



Product	Gastroscope	Colonoscope
Model	XZING-W200B	XZING-C200B
Field of view	140°	140°
Depth of field	3-100mm	3-100mm
Illumination	LED	LED
Working Length	1300mm	1300mm
Channel Inner Diameter	3.2mm	3.2mm
Insertion Tube Outer Diameter	11mm	14mm
Bending Angulation range	Up 180° Down 160° Left 160° Right 160°	Up 180° Down 180° Left 160° Right 160°
Direction from which accessories enters		
Direction from which the water jet enters		